# PROVOX<sup>®</sup> NiD<sup>™</sup>






**Rx**  
ONLY











**CE** 0413










42° C  
108° F  
RT  
2° C  
36° F

## **Prescription information**

**CAUTION:** United States Federal law restricts this device to sale, distribution and use by or on order of a physician or a licensed practitioner. The availability of this product without prescription outside the United States may vary from country to country.

## **Disclaimer**

Atos Medical offers no warranty - neither expressed nor implied - to the purchaser hereunder as to the lifetime of the product delivered, which may vary with individual use and biological conditions. Furthermore, Atos Medical offers no warranty of merchantability or fitness of the product for any particular purpose.

## **Patents and trademarks**

Provox<sup>®</sup> is a registered trademark owned by Atos Medical AB, Sweden. Provox<sup>®</sup> NiD<sup>™</sup> is a trademark of Atos Medical AB. For information about protective rights (e.g. patents), please refer to our web page [www.atosmedical.com/patents](http://www.atosmedical.com/patents).

Manufacturer; Hersteller; Fabrikant; Fabricant;  
 Produttore; Fabricante; Fabricante; Tillverkare; Pro-  
 ducent; Producent; Valmistaja; Framleiðandi; Tootja;  
 Ražotājs; Gamintojas; Výrobce; Gyártó; Výrobca;  
 Proizvajalec; Producent; Producător; Proizvođač;  
 Proizvođač; Κατασκευαστής; Производител;  
 Üretici; მწარმოებელი; Ἀρτοποιητής; İstehsalçı;  
 Производитель; Produzen; Pengilang; 製造業者;  
 제조사; 製造商; 制造商; היצרן; الجهة المصنعة

Date of manufacture; Herstellungsdatum; Productieda-  
 tum; Date de fabrication; Data di fabbricazione; Fecha  
 de fabricación; Data de fabrico; Tillverkningsdatum;  
 Fremstillingsdato; Produksjonsdato; Valmistuspäivä;  
 Framleiðsludagsetning; Valmistamise kuupäev;  
 Ražošanas datums; Pagaminimo data; Datum výroby;  
 Gyártás időpontja; Dátum výroby; Datum proizvodnje;  
 Data produkcji; Data fabricației; Datum proizvodnje;  
 Datum proizvodnje; Ημερομηνία κατασκευής; Дата на  
 производство; Üretim tarihi; წარმოების თარიღი;  
 Ἀρτοποιημάτων ἀμύσηρη; İstehsal tarixi; Дата  
 изготовления; Tanggal produksi; Tarikh pembuatan;  
 製造日; 제조일자; 製造日期; 制造日期;  
 תאריך ייצור; تاريخ التصنيع

Use-by date; Verfallsdatum; Houdbaarheidsdatum;  
 Date limite d'utilisation; Data di scadenza; Fecha  
 de caducidad; Data de validade; Används senast;  
 Anvendes inden; Brukes innen; Viimeinen käyttöpäivä;  
 Síðasti notkunardagur; Kölblik kuni; Derīguma ter-  
 miņš; Naudoti iki; Datum expirace; Felhasználhatóság  
 dátuma; Dátum najneskoršej spotreby; Rok uporab-  
 nosti; Data ważności; Data expirării; Rok uporabe;  
 Rok upotrebe; Ημερομηνία λήξης; Срок на годност;  
 Son kullanım tarihi; ვარგისობა; Ἐπισυνάληθη ἡμε-  
 ρῶν; Son istifadə tarixi; Срок годности; Gunakan  
 sebelum tanggal; Tarikh guna sebelum; 使用期限;  
 사용 기한; 使用截止日期; 使用截止日期;  
 תאריך انتهاء الصلاحية; تاريخ انتهاء الصلاحية

Batch code; Chargencode; Partijnummer; Code de lot;  
 Codice lotto; Código de lote; Código do lote; Batch-  
 kod; Batchkode; Batchkode; Eränumero; Lotunúmer;  
 Partii kood; Partijas kods; Partijos kodas; Kód šarže;  
 Sarzs kód; Kód šarže; Številka serije; Kod partii; Codul

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## Part I: General and user information

### Glossary

Clinician	Medical professional or licensed speech and language therapist / pathologist or clinical specialist nurse who is trained in voice rehabilitation procedures.
HME	Heat and Moisture Exchanger (Artificial Nose). Device that retains the heat and moisture in the exhaled air, which otherwise is lost when breathing through a tracheostoma.
PE-Segment	Pharyngo-esophageal segment. The part of the esophagus (food pipe) where sound is created by vibration of tissue when using a voice prosthesis.
Silicone	A latex-free material often used in medical devices
TE-puncture	Small artificial opening created between the trachea (wind pipe) and the esophagus (food pipe).
Tracheostoma	Breathing opening in the front of the neck, where the windpipe is connected to the skin (also called 'stoma').
Voice prosthesis	A one-way valve with retaining flanges inserted into a TE-puncture to allow speech by redirecting air to the esophagus (food pipe) while reducing the risk of food and liquids entering the trachea (wind pipe).

### Intended use

The Provox NiD voice rehabilitation system is intended for use in prosthetic voice rehabilitation after total laryngectomy only by patients who have been trained in the use of the device and, as assessed by the clinician who prescribes the device, have demonstrated the ability to understand and consistently follow the Instructions for Use without clinician supervision.

The Provox NiD is intended for single patient use.

### CONTRAINDICATIONS

- Provox NiD is NOT intended for immediate insertion at the time of tracheo-esophageal puncture.
- Provox NiD may be contraindicated for patients with bleeding disorders or who are undergoing anticoagulant treatment. These patient categories need to be carefully and regularly assessed by their clinician for the risk of bleeding from the puncture.
- A small tracheostoma may make it difficult or impossible to insert or replace a Provox NiD Voice Prosthesis through the stoma.

## Description of the device

### Provox NiD Voice Prosthesis and Inserter

The Provox NiD Voice Prosthesis is available in two outer shaft diameters (17 and 20 French) and several lengths.

Each voice prosthesis comes individually packaged with a non-sterile, reusable NiD Inserter.

#### Provox NiD Voice Prosthesis (Figure 1)

- (1) Silicone valve flap
- (2) Radio-opaque valve seat
- (3) Voice prosthesis shaft
- (4) Esophageal flange
- (5) Tracheal flange
- (6) Voice prosthesis tip
- (7) Safety strap
- (8) Safety medallion with thread
- (9) Size information (diameter and length between flanges). Available sizes: 6, 8, 10, 12, 14 and 18 mm.
- (10) Outer Diameter:  
Provox NiD 17: 5.67 mm (17 Fr)  
Provox NiD 20: 6.67 mm (20 Fr)

#### Provox NiD Inserter (Figure 2)

- (1) Inserter tip
- (2) Bendable inserter shaft
- (3) Safety strap retainer
- (4) Release lever
- (5) Handle

#### Materials:

Voice prosthesis: Silicone rubber and fluoroplastic.

Medallion and thread: Silicone rubber and polypropylene.

Inserter: Polypropylene.

### Accessories for use with Provox NiD Voice Prosthesis

The Provox NiD has a full range of accessories for the care and maintenance of the voice prosthesis. These accessories include: a Provox Dilator (Figure 3) for dilatation of the puncture prior to insertion; a Provox Brush (Figure 4) and a Provox Flush (Figure 5) for daily cleaning. The accessories are available in different diameters and sizes corresponding to the diameter and length of the voice prosthesis used.

### Recommended lubricants

The recommended water-soluble lubricants are:

KY Jelly®; SurgiLube®; Xylocaine®/(lidocaine) Jelly.

CAUTION: Lubricants that are not water-soluble should not be used at any time.

## WARNINGS

- **Never remove the safety strap.** Removal will severely increase the risk that your voice prosthesis dislocates and falls into the trachea.
- **Never remove the safety medallion without consulting your clinician first.** The medallion is important for your safety. It helps reduce the risk of injury in the event that your voice prosthesis dislocates and falls into the trachea.

- **Never remove the voice prosthesis from the puncture without immediately inserting a Provox Dilator or new voice prosthesis.** This may result in aspiration and/or shrinkage or closure of the puncture, which may necessitate dilatation or re-puncture by a clinician.

## **PRECAUTIONS**

### **To reduce infection risk:**

- Make sure that your hands are thoroughly clean and dry before using your hands in the area of your tracheostoma, and before cleaning your voice prosthesis and/or accessories.
- Keep all devices clean that may enter or make any contact with your tracheostoma and voice prosthesis.

### **To reduce risk of product damage:**

- Only use drinking water for cleaning. Never use detergents other than the ones your clinician has recommended. Detergents may weaken or damage the product or may be absorbed by the product which could cause tissue reactions.
- Only use recommended water-soluble lubricants. Oil-based lubricants (e.g., Vaseline) may weaken or damage the product and should never be used .
- Only use genuine Provox accessories for the handling and cleaning of your voice prosthesis.

### **To reduce risks when using the device in combination with other products:**

- If you use stoma attachment devices (e.g., adhesives, laryngectomy tubes, or stoma buttons), make sure you do not hook the tracheal voice prosthesis flange, safety strap, medallion, or medallion thread. This may lead to pulling the voice prosthesis out of the puncture.
- If you use a laryngectomy tube or stoma button, make sure that the device does not exert pressure on the voice prosthesis. This may lead to severe tissue damage and/or accidental swallowing of the voice prosthesis.

## **Adverse events**

### **Consult your clinician if:**

- Any changes occur in the appearance of the material of the prosthesis or the way it fits in the puncture.
- You frequently suffer from coughing during drinking, and cleaning the prosthesis does not resolve it. This may be an early sign of leakage through or around the prosthesis.
- Speaking becomes difficult (requires higher effort) and cleaning the prosthesis does not resolve it.
- There are any signs of inflammation or tissue change at or near the puncture tract. Symptoms may include, but are not limited to, pain, heat, redness, and swelling.
- You have bleeding or overgrowth of tissue around the device.
- You have persistent pain or discomfort in the region of the Provox NiD.
- You experience chronic cough, trouble breathing, or blood in the mucus. These can be signs of a serious health condition that requires medical attention.

If you have lost the device from the puncture, swallowed, or aspirated the voice prosthesis contact your clinician immediately.

# Instructions for use

## Inspecting the package and Provox NiD before use

Always carefully inspect the product prior to each use (i.e. before insertion), and replace if needed. Never use the product if the package is damaged or opened. Make sure that:

- the valve flap of the voice prosthesis closes against the valve seat (Fig. 1:1, 2). Improper closure can lead to aspiration of fluids, which may result in coughing and/or serious lung injury.
- the voice prosthesis is intact (i.e., no tears, cracks, or other structural damage) (Fig. 1:3-6).
- the safety strap is intact (Fig.1:7). Hold the tracheal flange of the voice prosthesis between two fingers and gently pull the safety strap (Figure 6). If the safety strap appears stiff, weak or brittle, or if it shows any kind of damage, do not use the voice prosthesis and obtain a replacement.
- the safety medallion and thread, if used, are intact and connected to the voice prosthesis (Fig. 1:8). Hold the tracheal flange of the voice prosthesis between two fingers and gently pull the medallion (Fig. 7). In addition, carefully inspect the openings where the thread is attached to the voice prosthesis and the medallion. If the safety medallion and/or thread are damaged, do not use the voice prosthesis.

## Using the Provox NiD to speak

When you close off your tracheostoma, you can direct the air from your lungs through the prosthesis into your esophagus. That air stream causes vibration of the tissues in your esophagus, which produces the sound of your voice. You can speak by using a Heat and Moisture Exchanger like the Provox HME, using a handsfree speaking device like Provox Freehands HME or by occluding your stoma directly with a finger.

An HME may improve pulmonary function and also facilitates stoma occlusion and speaking. A handsfree device enables speech without the need to close the tracheostoma by a finger.

## General preparation for maintenance

**Caution:** Before use or maintenance of the voice prosthesis or use of any of the accessories, make sure that you have read and understand the Instructions for Use accompanying the products. It is also recommended that you fully practice all maintenance procedures under the supervision of your clinician before you perform any voice prosthesis maintenance on your own.

Make sure that your hands are thoroughly clean and dry before using your hands in the area of your tracheostoma, and before cleaning your voice prosthesis and/or accessories.

Avoid swallowing while the puncture is open. This may lead to aspiration of saliva and/or coughing. To reduce the risk of such aspiration and/or coughing, try keeping your mouth open while you replace the voice prosthesis (this prevents you from swallowing).

**Do not maintain, remove or replace the voice prosthesis directly after a meal.** This may lead to gagging and/or vomiting.

It is recommended that you alternate between two voice prostheses in order to facilitate voice prosthesis maintenance and to reduce the risk of aspiration, shrinkage and/or closure of the puncture.

## Prepare

Have the following items ready: voice prosthesis of the same size and diameter, a Provox Dilator (if used), hypoallergenic tape (or Provox Silicone Glue), and a water-soluble lubricant. Also make sure that you have a magnifying mirror and sufficient light to illuminate the tracheostoma.

1. Attach an intact, clean Provox NiD onto the Provox NiD Inserter.
2. Secure the voice prosthesis. Attach the safety strap to the inserter by placing it along the inserter shaft into the strap retainer area.  
**Caution:** Close the strap retainer flap and make sure that it is locked (Fig. 11 and 12).
3. Put the voice prosthesis on the tip of the inserter. The safety strap should now be under slight tension (Fig. 13).
4. Lightly lubricate the tip of the voice prosthesis with water-soluble lubricant (see section “Recommended Lubricants”) (Fig. 14).

## Voice prosthesis removal

Prepare the Provox Dilator by lightly lubricating the tip with water-soluble lubricant (Fig. 8). Avoid touching any part of the Provox Dilator that will be inserted into the puncture.

To remove the voice prosthesis, first carefully loosen the safety strap from the skin while retaining the voice prosthesis in the puncture by gently pushing on it with a finger. Hold the safety strap of the voice prosthesis securely and pull the voice prosthesis carefully out of the puncture (Fig. 9). Never use the medallion or the medallion thread to pull out the voice prosthesis.

Gently insert the Provox Dilator into the puncture until you reach the stop in front of the Provox Dilator handle. Tape the medallion to the skin, away from the stoma (Fig. 10).

## Voice prosthesis insertion

If your puncture is situated at an angle within your stoma, you may bend the Provox NiD Inserter slightly (Fig. 15). Be aware that after bending the shaft you should not rotate the inserter when placing the voice prosthesis.

1. Carefully remove the Provox Dilator (if used) from the puncture (Fig. 16).
2. Hold the inserter firmly with one hand with the safety strap release lever pointing upwards.  
**Caution:** Avoid touching the release lever (Fig. 17).
3. Gently insert the Provox NiD Voice Prosthesis into the TE puncture, until the tracheal flange reaches the tracheal mucosa (Fig. 18, 19).
4. Gently move the inserter back and forth and if the inserter is not bent rotate it slightly (Fig. 20) to make sure that the esophageal flange has unfolded correctly. When the esophageal flange has unfolded, press the release lever to free the safety strap (Fig. 21).
5. Gently remove the inserter from the puncture while holding the safety strap of the voice prosthesis against the skin of your neck (Fig. 22). The voice prosthesis should stay properly in place.

6. Use a piece of tape or Provox Silicone Glue to attach the safety strap and safety medallion (if used) to the skin (Fig. 23).
7. Clean the Provox NiD Inserter and let dry (Fig. 24).

The insertion and/or removal procedures may cause temporary bleeding. Consult with your clinician immediately if such bleeding persists.

## Cleaning the Provox NiD

In order for you to speak with your voice prosthesis, the prosthesis needs to be clean so that air can pass through it. Regular cleaning of the prosthesis also helps to prolong device life.

### Cleaning the voice prosthesis after removal

Clean the used voice prosthesis using drinking water. Use the Provox Brush to remove all debris from both the inside and outside of the voice prosthesis. After cleaning, put the voice prosthesis on clean gauze and let it dry. Make sure that the voice prosthesis is stored at room temperature and that it is protected from heat and direct sunlight.

**CAUTION:** Avoid using detergents to clean the voice prosthesis other than those recommended by your clinician.

### Cleaning the voice prosthesis inside the puncture

**CAUTION: Never clean the voice prosthesis inside the puncture if the safety medallion has been removed without a clinician's recommendation, or if the medallion is damaged.** The medallion is important to your safety. It helps reduce the risk of injury in the event that your voice prosthesis dislocates and falls into the trachea. **The strap must be secured to the skin when cleaning the VP inside the puncture.**

Before cleaning the voice prosthesis, make sure the safety strap of the voice prosthesis is secured to the skin. Additionally, when cleaning the voice prosthesis inside the puncture, secure the safety strap by pressing it against the skin with a finger.

Prior to use, carefully read the Instructions for Use that accompanies Provox Brush and Provox Flush.

- Brushing the prostheses with the Provox Brush twice a day and after each meal removes mucus and food remnants from the prosthesis.
- Flushing the prosthesis with the Provox Flush also helps to clean debris and mucus from the prosthesis, which can help increase the life of the device.

**Note:** The Provox Flush is intended for use only by patients who have been assessed by the clinician to be able to consistently follow the instructions for use without clinician supervision.

- Some dietary measures, like the daily intake of yogurt or butter milk containing lactobacilli, are considered to be helpful against excessive *Candida* growth.

**CAUTION:** Only use genuine Provox accessories that are intended for use with Provox NiD when cleaning your prosthesis. Other devices may cause personal injury or damage to the voice prosthesis.

## Device lifetime

If the voice prosthesis functions properly and has no structural damage, it may be used for a maximum period of 6 months. If the voice prosthesis does not function properly, and cleaning does not solve the problem, discard the device and replace it with a new one immediately.

## Disposal

Always follow medical practice and national requirements regarding biohazard when disposing of a used medical device.

## Troubleshooting: - Events by symptom

### **Coughing (triggered by drinking or eating) while the prosthesis is in place.**

**Most common reasons:** Leakage through or around the prosthesis.

**What to do:** Clean the prosthesis with the Provox Brush and/or Flush. There may be food remnants preventing the valve from closing. If cleaning does not solve the problem, discard the device and replace it with a new one immediately. If replacement of the voice prosthesis does not improve symptoms, you may have leakage around the voice prosthesis. Consult with your clinician immediately.

### **Speaking becomes more and more difficult.**

**Most common reasons:** Valve obstruction, tissue swelling in the throat (e.g. during radiotherapy).

**What to do:** Clean the prosthesis with the Provox Brush and/or Flush. There may be food remnants increasing the airflow resistance. If cleaning does not improve symptoms, stop using the voice prosthesis and replace it with a new one immediately. If replacement of the voice prosthesis does not improve symptoms, consult with your clinician immediately.

### **Insertion of the voice prosthesis into the puncture is difficult.**

**Possible causes:** Shrinkage of puncture. Swelling of tissue, excessive tissue, or deformation of the puncture. Not enough lubricant on the voice prosthesis.

**What to do:** Stop insertion immediately, and instead gently insert a Provox Dilator of the correct size. With the dilator in place, wait for a few minutes and then retry inserting the voice prosthesis with some lubricant on its tip. If the difficulties continue, leave the dilator in place and consult with your clinician immediately.

**CAUTION:** DO NOT leave the puncture empty after removal of your voice prosthesis; always insert the dilator or a voice prosthesis. The puncture may shrink, requiring dilatation or re-puncture by a clinician.

### **Inability to insert dilator.**

**Possible causes:** Shrinkage of puncture. Deformation of the puncture:

**NEVER** use force to insert the dilator. Consult with your clinician immediately.

### **After insertion, voice prosthesis slides out from the puncture.**

**Possible causes:** Too much lubricant used during voice prosthesis insertion. Enlarged puncture.

**What to do:** If too much lubricant was used, wipe off the excess lubricant and reinsert the voice prosthesis. Secure the safety strap to the skin and hold the voice prosthesis

with a finger in place. A warm drink may help to dissolve the lubricant on the esophageal side.

If the puncture is enlarged, you may also notice leakage, even if the voice prosthesis is clean and intact. This can be observed as leakage around the voice prosthesis, instead of through the voice prosthesis. Consult with your clinician immediately if this is the case.

**Pain, heat, swelling, and/or rash in the area of the puncture or stoma.**

**Possible causes:** Tissue inflammation and/or infection.

**What to do:** Consult your clinician immediately. If you experience persistent bleeding from the puncture, seek immediate medical attention.

**Voice prosthesis cannot be seen from outside; safety medallion thread leads into the puncture.**

**Possible causes:** Accidental dislodgement of the voice prosthesis in the direction of the esophagus.

**What to do:** Grasp the safety strap and pull the voice prosthesis carefully back in place. If the strap is not visible, grasp the thread of the safety medallion as close as possible to the puncture. Gently pull the thread in the direction of the puncture. As soon as the safety strap or voice prosthesis flange appears, grasp it and pull the voice prosthesis correctly in place.

If pulling the thread is not successful, gently dilate the puncture by inserting the Provox Dilator (while holding the thread that remains in the puncture), remove the Dilator and retry the pulling procedure described above. If still not successful, gently insert the Provox Dilator again (while holding the thread that remains in the puncture) to avoid puncture shrinkage. Consult with your clinician immediately.

**Voice prosthesis falls into trachea and safety medallion is at the level of the tracheostoma:**

**What to do:** Breathe calmly and remove the voice prosthesis from the trachea by gently pulling the medallion up. Insert a Provox Dilator into the puncture, and prepare and insert a clean voice prosthesis. If the voice prosthesis falls into the trachea more than once, consult with your clinician immediately.

**Voice prosthesis is missing and safety medallion is removed:**

**Possible causes:** Accidental swallowing or accidental aspiration of the voice prosthesis. Voice prosthesis may have been dislodged and misplaced outside the body or voice prosthesis is hidden under tissue in the puncture tract.

**What to do:** If the voice prosthesis cannot be coughed up at once or is not found, gently try to insert a Provox Dilator. Breathe calmly and seek immediate medical attention to find and if necessary remove the voice prosthesis. The blue ring inside the voice prosthesis can be detected by X-ray.

## **User assistance information**

For additional help or information please see back cover of the Instructions for Use for contact information.

## **Part II: Prescriber information**

### **User selection**

#### **General requirements**

Potential users of the Provox NiD must have sufficient cognitive ability, manual dexterity, and visual acuity to understand and follow both the Instructions for Use supplied

with this device and the directions furnished by the clinician in the exercise of his or her independent medical judgment. It is the responsibility of the prescribing clinician to assess and ensure complete user compliance with all warnings, precautions and instructions appropriate to the use of this device. Also consider aging, it is recommended that the prescribing clinician determines each patient's ability and understanding regarding use of the product at least once a year. If the patient no longer fulfills the requirements consider conversion to an indwelling voice prosthesis.

If you doubt in any way that your patient fulfills all of the above criteria, including without limitation that he or she has fully understood the information you have provided in the exercise of your independent medical judgment, do not prescribe the Provox NiD voice rehabilitation system.

### **Removal of safety medallion**

If the patient wishes to use the voice prosthesis without the safety medallion it is the responsibility of the clinician to determine if the retention of the voice prosthesis in the TE puncture is sufficient to do so. The retention of the voice prosthesis must also be checked regularly at an interval determined by the clinician to ensure that changes in the puncture have not compromised the retention of the voice prosthesis.

### **User training and counseling**

Each user must be fully instructed as to the procedures, risks, and other information deemed necessary by the clinician, including adverse events, their symptoms, and all necessary and appropriate actions to be taken by the user in the event that such adverse events occur.

## **Choosing the right voice prosthesis**

### **Diameter**

Provox NiD prostheses are available in outer shaft diameters of 17 and 20 Fr (5.67 and 6.67 mm, respectively). Although both voice prosthesis types are optimized for minimum airflow resistance at normal speech loudness, the 20 Fr diameter voice prosthesis has lower airflow resistance, which may allow for easier voicing.

Use a Provox NiD Dilator of the size corresponding to the selected voice prosthesis to ensure that the pathway of the voice prosthesis is the appropriate size prior to insertion.

To allow for an increase from 17 Fr to 20 Fr, dilate the puncture first using the Provox NiD 20 Dilator.

**Note:** The Provox Dilators will dilate the puncture tract slightly more than the corresponding voice prosthesis size, i.e., the Provox Dilator 17 has a 19 Fr shaft diameter, and the Provox Dilator 20 has a 22 Fr shaft diameter, in order to facilitate voice prosthesis insertion.

### **Length**

A too short or long voice prosthesis may cause severe tissue damage and/or problems with speech.

Always check the appropriate length of the voice prosthesis. Use the old voice prosthesis in place to check the size required or use a Provox Measure to determine the appropriate size.

## **Additional information**

### **Compatibility with antifungal medicine**

Chemical substances may influence the material properties of the device or may be absorbed by the material and subsequently released at a later time. Therefore, the introduction of antifungal medicines or drugs directly to or within close proximity of the voice prosthesis should be carefully assessed.

Laboratory tests show only minor negative influence on the function of the Provox NiD Voice Prosthesis and components when using the following antifungal medications:

- Nystatin (Mycostatin<sup>®</sup>, Bristol-Myers Squibb, NY, USA).
- Fluconazole (Diflucan<sup>®</sup>, Pfizer, NY, USA).

### **Compatibility with cleaning detergents**

The preferred method to clean the Provox NiD is to use drinking water without any detergents, as described in the User Manual, chapter “Cleaning the Provox NiD”.

However, laboratory tests show only moderate negative influence on the function of the Provox NiD Voice Prosthesis and components with repeated short-term exposure no more than 15 minutes per day to: (listed in the order of increasing negative effect on device function).

- Efferdent Anti Bacterial Denture Cleanser (Pfizer Consumer Healthcare, NJ, USA),
- Ethanol 70%,
- Hydrogen peroxide (3%) solution.

## **Reporting**

Please note that any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the national authority of the country in which the user and/or patient resides.

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