FREEVENT

Dressing Softfoam


Figure 1


Figure 2

Disclaimer

Atos Medical offers no warranty—neither expressed nor implied—to the purchaser hereunder as to the lifetime of the product delivered, which may vary with individual use and biological conditions. Furthermore, Atos Medical offers no warranty of merchantability or fitness of the product for any particular purpose.

	Length (mm)	Width (mm)	Thickness (mm)
Freevent Dressing Softfoam L	80	80	5
Freevent Dressing Softfoam S	50	80	5
Freevent Dressing Softfoam Slim L	80	80	2,5

Contents

EN - ENGLISH.....	5
DE - DEUTSCH.....	7
NL - NEDERLANDS.....	8
FR - FRANÇAIS.....	10
IT - ITALIANO	11
ES - ESPAÑOL.....	13
PT - PORTUGUÊS.....	14
SV - SVENSKA	16
DA - DANSK.....	17
NO - NORSK.....	18
FI - SUOMI.....	20
LT - LIETUVIŲ K.	21
CS - ČESKY	22
HU - MAGYAR	24
SK - SLOVENČINA.....	25
PL - POLSKI	27
RO - ROMÂNĂ	28
HR - HRVATSKI.....	30
SR - SRPSKI	31
TR - TÜRKÇE.....	33
RU - РУССКИЙ.....	34
AR - العربية.....	37

1. Intended Use

Freevent Dressings are sterile single-use dressings that absorb secretion and provide protection for the skin around the tracheostomy tube for patients with a tracheostoma.

2. Warnings

- When changing the dressing, be cautious in order to prevent accidental de-cannulation.

3. Precautions

- Freevent Dressings may be used on irritated skin. However, if skin irritation persists for more than 14 days, stop using the dressing and consult your clinician.
- The dressing should be changed at least every 24 h, or more often if needed, to prevent skin irritation.
- Do not use the dressing if the pouch is damaged.
- Do not modify the dressing by e.g. cutting.
- Applying the wrong side of the dressing towards the skin may cause skin irritation.
- The dressings may irritate the skin. Stop using the dressing if skin irritation develops and consult your clinician.

4. Description of the Device

Freevent Dressing Softfoam is made of soft polyurethane foam that provides medium to high absorption of secretion.

5. How to Use

5.1. Operating

Freevent Dressing Softfoam has a pre-perforated slit that can be torn open (see Figure 1).

1. Apply the dressing with the foam side against the skin and shiny side up. The dressing should be placed between the tracheostomy tube and the skin with the slit facing upwards (see Figure 2).
2. When using a dressing without opening the slit, insert the tracheostomy tubes through the opening of the dressing.
3. Insert the tracheostomy tube in the tracheostoma and fix it with a tube holder.

6. Disposal

Freevent Dressings shall be changed at least every 24 h. Always follow medical practice and national requirements regarding biohazard when disposing of a used medical device.

7. Reporting

Please note that any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the national authority of the country in which the user and/or patient resides.

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Atos and the Atos logo are trademarks of Coloplast A/S.
© 2024 Coloplast A/S. All rights reserved.
Atos Medical AB